



## BUCCUTITE™ CROSSLINKING PRINCIPLES

### Conjugate Proteins In 2 Simple Mixing Steps.


\*NO catalyst required



1. Desired antibody is activated with MTA.



2. Activated IgG-Buccutite™ MTA and preactivated HRP-Buccutite™ FOL are mixed



3. HRP-labeled antibody is ready to use for desired immunoassay application.

\*Buccutite™ Conjugation using Buccutite™ Peroxidase (HRP) Antibody Conjugation Kit (Catalog: 5503)


## PRODUCT ORDERING INFORMATION FOR BUCCUTITE™ LABELLING KITS

Cat #	Product Name	Abs/Em (nm)
5503	Buccutite™ Peroxidase (HRP) Antibody Conjugation Kit *Optimized for Labeling 100 µg Protein*	N/A
1313	Buccutite™ Rapid APC Antibody Labeling Kit *Microscale Optimized for Labeling 2w5 µg Antibody Per Reaction*	651/662
1311	Buccutite™ Rapid APC Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	651/662
1322	Buccutite™ Rapid APC-Cy5 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	651/670
1350	Buccutite™ Rapid APC-Cy5.5 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	651/700
1320	Buccutite™ Rapid APC-Cy5.5 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	651/700
1351	Buccutite™ Rapid APC-Cy7 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	651/780
1321	Buccutite™ Rapid APC-Cy7 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	651/780
1347	Buccutite™ Rapid APC-iFluor™ 700 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	651/713
1319	Buccutite™ Rapid APC-iFluor™ 700 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	651/713
1312	Buccutite™ Rapid PE Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	565/575
1310	Buccutite™ Rapid PE Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	565/575
1340	Buccutite™ Rapid PE-Cy5 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	565/670
1341	Buccutite™ Rapid PE-Cy5.5 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	565/700
1316	Buccutite™ Rapid PE-Cy5.5 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	565/700
1342	Buccutite™ Rapid PE-Cy7 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	565/700
1317	Buccutite™ Rapid PE-Cy7 Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	565/780
1353	Buccutite™ Rapid PerCP Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	482/677
1325	Buccutite™ Rapid PerCP Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	482/677
1343	Buccutite™ Rapid PE-Texas Red Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 25 µg Antibody Per Reaction*	565/600
1318	Buccutite™ Rapid PE-Texas Red Tandem Antibody Labeling Kit *Microscale Optimized for Labeling 100 µg Antibody Per Reaction*	565/600
1315	Buccutite™ Rapid Protein Crosslinking Kit *Microscale Optimized for Crosslinking 100 µg Antibody Per Reaction*	N/A

For more information contact us at 1-(800)-990-8053 or visit our website at [www.aatbio.com](http://www.aatbio.com)

# BUCCUTITE™ BIOCONJUGATION TECHNOLOGY

Rapid Protein Crosslinking Method for Labeling and Modifying Antibodies



 **AAT Bioquest®**

## CROSSLINKING METHODS

Crosslinking is a technique that chemically joins two or more proteins together, such as an antibody and an enzyme, by a covalent bond. Crosslinking reagents consist of reactive moieties that modify and attach to specific functional groups on proteins. A common approach for conjugating two biomolecules utilizes a small, heterobifunctional crosslinker called succinimidyl 4-(N-maleimidomethyl) cyclohexane-1-carboxylate (SMCC). SMCC crosslinkers contain amine-reactive NHS ester moieties and thiol-reactive maleimide moieties to crosslink proteins. However, a drawback to this approach is the many experimental conditions that must be monitored to perform a successful conjugation. Difficulties associated with SMCC conjugations are maintenance of a narrow pH range, self-polymerization, and hydrolytic degradation of reactive moieties. To address these concerns, AAT Bioquest has developed an optimized and robust Buccutite™ crosslinking technique. Buccutite™ crosslinking technology provides a simplistic and efficient approach for conjugating proteins. It can be completed in half the time and with less stringent parameters yielding highly stable and easy to handle conjugated proteins.

## BUCCUTITE™ vs SMCC

### TOTAL OPERATION TIME:

- **BUCCUTITE™ CHEMISTRY:** 2 HOURS
- **SMCC CHEMISTRY:** 4-8 HOURS

### MINIMUM SAMPLE CONCENTRATION:

- **BUCCUTITE™ CHEMISTRY:** ≥ 0.5 MG/ML
- **SMCC CHEMISTRY:** 0.5 - 5.0 MG/ML

### OPTIMAL CONJUGATION pH:

- **BUCCUTITE™ CHEMISTRY:** 5 - 9
- **SMCC CHEMISTRY:** 7.0 - 7.5

### CONJUGATION YIELD:

- **BUCCUTITE™ CHEMISTRY:** 50-60 %
- **SMCC CHEMISTRY:** ≤ 30 %

## ADVANTAGES

### FAST

Buccutite™ conjugation can be completed in 1-2 hours under extremely mild conjugation conditions. No purification step required with 25µg kits.

### ROBUST


Buccutite™ conjugation can be run in a broad pH range of 5-9 and temperature.

### HIGHLY STABLE

Buccutite™ linker-activated macromolecules are very stable, can be stored at 4 °C for more than 24 months.


### HIGH YIELD

Buccutite™ conjugation gives much higher yield than other existing methods under the same conditions. Eliminating homo-crosslinking of proteins.



Tubulins were imaged with RPE- goat-anti-mouse IgG conjugate in HeLa cells. Tubulins were stained using mouse anti-α-tubulin antibody, and visualized with red fluorescent RPE-goat-anti-mouse IgG conjugate prepared by Buccutite™ crosslinking technology as described above.

## RESULTS



Flow cytometry analysis of HL-60 cells stained with 1µg/ml Mouse IgG control (Green) or with 1µg/ml mouse Anti-Human HLA-ABC (W6/32 mAb) (Red) and then followed by Goat Anti-Mouse IgG-RPE conjugate prepared with Buccutite™ Rapid RPE Antibody Labeling Kit (Cat#1310). The fluorescence signal was monitored using ACEA NovoCyte flow cytometer in the RPE channel.

## AAT BIOQUEST'S CUSTOM BIOCONJUGATION SERVICES

**\*TOGETHER WE SHINE<sup>SM</sup>\***

AAT Bioquest offers same-day custom conjugation of proteins or antibodies with a wide array of labels, such as biotin, HRP and over 20 fluorophores. Our services are designed to produce high-quality results with a fast turnaround. We take pride in supplying scientists with flexible tools that help expand their research needs.

### Features of AAT Bioquest's custom bioconjugation services include:

- Guaranteed quality with 95% purity rating
- Affordable price tiers that work with your budget
- Scalable service with minimum order of 50 µg
- Same day order fulfilment means your conjugate ships the day you order it

**\*NOTE: APC and PE custom conjugation services require 3 working days**

You may supply your own protein or antibody or chose from 3000+ monoclonal and polyclonal antibodies in our catalog. Contact us for a quote today, and let us put our expertise to work for you!