


Agility[®]

*Fully automated ELISA processor
compatible with TestLine SmartKits[®]*





The instrument redefining ELISA processing with major impact on automation, ease of use, reagents loading, hands on time, productivity and intelligence. Now compatible also with TestLine SmartKits® for the most comprehensive ELISA utilization.

Agility® significantly reduces the total assay time as well as loading of reagents. The system allows up to 16 SmartKits® to run simultaneously as well as to keep 12 plates on board at one time. All consumables for each assay are already prepackaged in the TestLine SmartKit® for the maximum user comfort and safety.

The new dimension in ELISA processing
The true solution for walkaway performance

Automation

- ▶ Up to 16 SmartKits® carriers stored on-board for simultaneous runs
- ▶ Flexible throughput allows up to 12 plates on-board at once
- ▶ Up to 200 sample tubes on-board with additional continuous load capability
- ▶ Utilizes three precision robotic arms to obtain maximum process efficiency
- ▶ Comprehensive monitoring system measuring consumables levels



Ease-of-use & hands on time

- ▶ Intuitive, easy-to-use interface
- ▶ Eliminates most of ELISA's labor-intensive, front-end setup
- ▶ Reduces hands-on time by two-thirds of typical open systems
- ▶ Barcoding eliminates manual data input
- ▶ Comprehensive LIS integration
- ▶ Continuous sample loading allows operator to begin loading and running microplates as they are ready, instead of all at once



Agility™ Hands on Time

12:55

Minutes



Traditional Open Systems
Hands on Time

36:45

Minutes



Reagent loading

- ▶ Eliminates nearly all manual liquid-transfer steps
- ▶ Free up significant labor time and allows for multitasking
- ▶ Illuminated, color-coded access positions allow simple loading of all kits, tips, microplates and other consumable items, preventing loading errors
- ▶ 2D barcodes printed on kits provide accurate and secure information
- ▶ Inventory tracking allows for maximum materials efficiency

Safety

- ▶ Pre-packed TestLine SmartKits® with printed barcodes
- ▶ Sample tube barcode scan
- ▶ Automatic tube type recognition
- ▶ Automatic strip check
- ▶ Color-coded access positions

Agility® Parameters

Dimensions

Width:	<1250 mm
Depth:	<900 mm
Height:	<1230 mm
Footprint:	<1200 x 650 mm
Bench Weight:	213 kg (max.)
Ship Weight:	296 kg (max.)
Noise:	Noise output <80 dB

Power Requirements

Voltage:	100 – 240V automatic switching
Frequency:	50/60 Hz
Power Consumption:	Typically 400 VA (online UPS required)

Reader Specifications

Photometric Range:	0.000 to 3.500 OD
Spectral Range:	405 nm to 690 nm
	±0.010 OD at 0.000 to 0.500 OD
Precision:	<1% CV at 0.501 to 2.000 OD
	<1.5% CV at 2.001 to 2.500 OD
Accuracy:	±0.01 or 2.5%
	(0.000 to 3.000 OD) whichever is greater
Read Time:	<30 seconds, single wavelength
	<50 seconds, dual wavelength

Washer Specifications

Manifold Configuration:	8-way wash head
Programmable Volumes:	50-999 µl
Wash Containers:	4 wash bottles at 3.0 L, with quantitative level sensing
Clean Container:	1 wash bottle at 3.0 L, with quantitative level sensing
Waste Container:	10 L with quantitative level sensing
Residual Wash Volume:	<3 µl per well with dual-axis sweep in a at bottom well
Dispense Precision:	≤5% CV (with 300 µl in a 96 well plate)

Incubator Specifications

Number of Incubators:	Up to 12 (6 Ambient / 6 Elevated)
Temperature Range:	RT + 4° C to 45° C (Elevated Incubator) RT + 4° C (Ambient Incubator)

Temperature Accuracy:	± 1° C
Shaking:	14 Hz periodic or continuous

System Specifications

Number of Plates:	12 (additional continuous load capability)
Number of Sample Tubes:	200 (additional continuous load capability)
Number of Reagents:	Up to 16 different SmartKit™ Reagent Packs simultaneously (additional packs can be added using continuous load capability)
Number of Pipettes:	2 (sample + reagent)

Sample Pipetting Specifications

Sample Tip Size:	300 µl
Sample Pipetting Volume:	10-300 µl
Time to Dispense:	11 minutes (typical) 50 µl of 96 samples to plate from sample tubes or deep well plates
Sample Pipetting Precision:	≤3% CV (at any operating volume above 10 µl)
Sample Pipetting Accuracy:	± 2% of target volume
Dilution Range:	1 part in 199 one-stage dilution, 1 part in 39,601 two-stage dilution
Number of Sample Tips:	20 racks of 112 tips (additional racks can be added using continuous load capability)
Sample Tube Dimensions:	10-17 mm diameter external dimension, 45-100 mm depth

Reagent Pipetting Specifications

Reagent Tip Size:	1200 µl
Reagent Pipetting Volume:	20-1200 µl
Number of Reagent Tips:	1 rack of 98 tips (additional racks can be added using continuous load capability)
Reagent Pipetting Precision:	≤3% CV at 10 shots at any volume in operating range
Reagent Pipetting Accuracy:	±2% of target volume (single shot mode)

Regulatory Compliance

Certified:	ISO13485
Electromagnetic Perturbation:	EMC Directive 2004/108/EC; EN 61326-2-6:2006, IEC 61326-1:2006
Electrical Safety:	BS EN 61010-1:2001, IEC 61010-1:2001, IEC 61010-2-101, ETL marked
Lead-free:	All components are RoHS compliant
IVDD:	CE marked per IVD 98/79/EC Directive

TestLine Clinical Diagnostics s.r.o.
Krizikova 68
612 00 Brno
Czech republic

ImTec
DIAGNOSTICS NV

