

Overview of Rapid, CLIA, RT-PCR & ELISA tests

Rapid Antigen Tests


| Catalog No. | Description | Sample Type | Target | Size per package | Information |
|-------------|--|---|--|------------------|---|
| RT2952 | AMP SARS-CoV-2 Ag | Naso- or Oropharyngeal or anterior nasal swab | Nucleocapsid protein antigen | 25 cassettes | AMP Rapid Test SARS-CoV-2 Ag |
| RT3952 | AMP SARS-CoV-2 Ag | Sputum | Nucleocapsid protein antigen | 20 cassettes | AMP Rapid Test Sputum |
| RT2962 | AMP SARS-CoV-2 Ag and Influenza A+B | Nasopharyngeal swab | Nucleocapsid protein antigen, Influenza A and/or B | 25 cassettes | IFU SARS-CoV-2 Ag & A+B |
| STCOV-101 | COVID-19 Saliva test | Saliva | Specific antibody-antigen reaction | 1 cassette | IFU COVID-19 Saliva test. Noviral |
| RADCOV7-20T | COVID-19 Antigen Rapid test | Saliva | SARS-CoV-2 antigens | 20 cassettes | COVID-19 Antigen Rapid test. Prestige |
| 500220 | SARS-CoV-2 & Influenza A+B Combo Antigen | Nasal/Oropharyngeal swab | Nucleocapsid protein antigen | 20 cassettes | IFU StrongStep® Combo antigen |

Rapid Antibody Tests

| Catalog No. | Description | Sample Type | Target | Size per package | Information |
|--------------|-----------------------------|------------------------------|---------------------------------|------------------|---|
| ABT-IDT-B352 | COVID-19 IgG/IgM Cassette | Whole blood, Serum or Plasma | Combination of COVID-19 antigen | 50 cassettes | IFU COVID-19 IgG/IgM Cassette Accubiotech |
| INGM-MC42 | COVID-19 IgG/IgM Rapid Test | Whole blood, Serum or Plasma | Combination of COVID-19 antigen | 25 cassettes | IFU COVID-19 IgG/IgM Rapid Test Cassette Surescreen |

CLIA analyzer kits

| Catalog No. | Description | Target protein | Size | Information |
|-------------|---------------------------|---|---------------|--|
| 130219016M | MAGLUMI™ 2019-nCoV IgM | Spike & nucleocapsid protein | 100 tests/kit | IFU MAGLUMI™ 2019-nCoV IgM (CLIA) |
| 130219015M | MAGLUMI™ 2019-nCoV IgG | Spike & nucleocapsid protein | 100 tests/kit | IFU MAGLUMI™ 2019-nCoV IgM (CLIA) |
| 130219026M | MAGLUMI™ SARS-CoV-2 Ag | Nucleocapsid protein | 100 tests/kit | IFU MAGLUMI™ SARS-CoV-2 Ag IgG (CLIA) |
| 130219017M | MAGLUMI™ SARS-CoV-2 S-RBD | Spike protein receptor-binding domain (RBD) | 100 tests/kit | IFU MAGLUMI™ SARS-CoV-2 S-RBD IgG (CLIA) |



ELISA Kits


| Catalog No. | Description | Target protein | Size |
|-------------|---------------------------|--|----------------|
| CoNA96 | EIA COVID-19 NP IgA | Nucleocapsid recombinant antigen | 96 wells plate |
| CoNG96 | EIA COVID-19 NP IgG | Nucleocapsid recombinant antigen | 96 wells plate |
| CoNM96 | EIA COVID-19 NP IgM | Nucleocapsid recombinant antigen | 96 wells plate |
| SK-CoNA96* | SmartEIA COVID-19 NP IgA | Nucleocapsid recombinant antigen | 96 wells plate |
| SK-CoNG96* | SmartEIA COVID-19 NP IgG | Nucleocapsid recombinant antigen | 96 wells plate |
| SK-CoNM96* | SmartEIA COVID-19 NP IgM | Nucleocapsid recombinant antigen | 96 wells plate |
| CoRA96 | EIA COVID-19 RBD IgA | Recombinant antigen Receptor-binding domain, a subunit of the Spike S1 protein | 96 wells plate |
| CoRG96 | EIA COVID-19 RBD IgG | Recombinant antigen Receptor-binding domain, a subunit of the Spike S1 protein | 96 wells plate |
| CoRM96 | EIA COVID-19 RBD IgM | Recombinant antigen Receptor-binding domain, a subunit of the Spike S1 protein | 96 wells plate |
| SK-CoRA96* | SmartEIA COVID-19 RBD IgA | Recombinant antigen Receptor-binding domain, a subunit of the Spike S1 protein | 96 wells plate |
| SK-CoRG96* | SmartEIA COVID-19 RBD IgG | Recombinant antigen Receptor-binding domain, a subunit of the Spike S1 protein | 96 wells plate |
| SK-CoRM96* | SmartEIA COVID-19 RBD IgM | Recombinant antigen Receptor-binding domain, a subunit of the Spike S1 protein | 96 wells plate |
| COVG0940 | SARS-CoV-2 (COVID-19) IgG | Nucleocapsid recombinant antigen | 96 wells plate |
| COVA0940 | SARS-CoV-2 (COVID-19) IgA | Nucleocapsid recombinant antigen | 96 wells plate |
| COVM0940 | SARS-CoV-2 (COVID-19) IgM | Nucleocapsid recombinant antigen | 96 wells plate |

*SmartEIA kits are designed for automated processing using the Agility® analyser.

Microblot-Array Kits

| Catalog No. | Description | Antigens | Size |
|-------------|------------------------------|--|----------|
| CoVAMA96* | Microblot-Array COVID-19 IgA | Nucleocapsid, RBD, Spike S2, Envelope protein, ACE2, PLPro protein, MERS-CoV, SARS-CoV, HCoV229E NP and HVoV NL63 NP | 96 tests |
| CoVGMA96* | Microblot-Array COVID-19 IgG | Nucleocapsid, RBD, Spike S2, Envelope protein, ACE2, PLPro protein, MERS-CoV, SARS-CoV, HCoV229E NP and HVoV NL63 NP | 96 tests |
| CoVMMA96* | Microblot-Array COVID-19 IgM | Nucleocapsid, RBD, Spike S2, Envelope protein, ACE2, PLPro protein, MERS-CoV, SARS-CoV, HCoV229E NP and HVoV NL63 NP | 96 tests |

*Microblot-Array is performed as an ELISA but read as a Blot



RT-PCR Kits

| Catalog No. | Description | Target gene | Rxn |
|----------------------|---|---|----------------|
| PCOV6050TC | NovaType SARS-CoV-2 Detect & ID | S gene | 288 or 96 |
| DHUV02610 | ViroReal® Kit RT-LAMP SARS-CoV-2 | ORF1ab region | 100 or 500 |
| DHUV02413 | ViroReal® Kit SARS-CoV-2 Multiplex | N & RdRp genes | 100 or 500 |
| DHUV02513 | ViroReal® Kit SARS-CoV-2 and Influenza A+B | N gene, matrix protein gene, hemagglutinin gene | 100 or 500 |
| DHUV02353/ DHUV02313 | ViroReal® Kit SARS-CoV-2 & SARS | N gene | 50 or 100 |
| 3239 | gb SARS-CoV-2 Variant SA/BR, UK | S & RdRP genes | 100 or 500 |
| 3231 | gb SARS-CoV-2 Multiplex | E & RdRP genes | 50, 200 or 500 |
| 3234 | gb SARS-CoV-2 Multiplex EndoC | E & RdRP genes | 100 or 500 |
| 3232 | gb SARS-CoV-2 Combi | E & RdRP genes | 100 or 500 |
| 3235 | gb SARS-CoV-2 Combi EndoC | E & RdRP genes | 100 or 500 |
| 3233 | gb SARS-CoV-2 Influenza A/B | E & RdRP genes, gene M | 100 or 500 |
| 3236 | gb SARS-CoV-2 Influenza A/B EndoC Coming soon! | | |

RT-PCR Duplex Kits (generi biotech)

| Catalog No. | Description | Target gene | Rxn | IFU |
|-------------|--|-------------|------------|--|
| 3229 | gb Sarbeco N (primary test) | N gene | 100 or 500 | IFU gb Sarbeco N (primary test) |
| 3227 | gb Sarbeco E (primary test) | E gene | 100 or 500 | IFU gb Sarbeco E (primary test) |
| 3228 | gb SARS-CoV-2 RdRP (confirmation test) | RdRP gene | 100 or 500 | IFU gb SARS-CoV-2 RdRP (confirmation test) |
| 3230 | gb SARS-CoV-2 N (confirmation test) | N gene | 100 or 500 | IFU gb SARS-CoV-2 N (confirmation test) |