

Backbone panels

Panel	Reactiviteit	Targets
Basic B Lymphocyte Backbone Panel	Human	CD4 PerCP-Cy5.5 CD19 PE CD20 redFluor 710 CD27 PE-Cy7 CD38 APC CD45 violetFluor 450
Basic T Lymphocyte Backbone Panel	Human	CD3 FITC CD4 violetFluor 450 CD8a PE CD45 PerCP-Cy5.5
Basic T Lymphocyte Backbone Panel	Mouse	CD3 violetFluor 450 CD4 APC-Cy7 CD8a redFluor 710 CD25 PE-Cy7 CD44 PerCP-Cy5.5 CD62L APC
Dendritic Cell Backbone Panel	Mouse	CD8a PE-Cy7 CD11b APC-Cy7 CD11c PerCP-Cy5.5 CD45R APC Ly-6G PE MHC Class II FITC
Myeloid Precursor Backbone Panel	Mouse	CD3 redFluor 710 CD11b violetFluor 450 CD11c PerCP-Cy5.5 CD19 redFluor 710 CD45R redFluor 710 LY-6G PE-Cy7 TER-119 FITC
T Lymphocyte Subset Backbone Panel	Mouse	CD4 PE CD3e PE-Cy7 CD8a APC
T Lymphocyte Subset Backbone Panel	Mouse	CD3e PE-Cy7 CD4 PE CD8a FITC
DryFlowEx ASC Screen Kit	Human	CD45 Pacific Blue IgD FITC CD27 PE CD24 PerCP-Cy5.5 CD19 PE-Cy7 CD21 APC CD38 APC-Cy7
DryFlowEx ACT T Screen Kit	Human	CD45 Pacific Blue CD8 Pacific Orange CD4 FITC PD-1 PE HLA-DR PerCP-Cy5.5 CD3 PE-Cy7 CXCR5 APC CD38 APC-Cy7

redFluor and violetFluor are trademarks of Tonbo Biosciences, Cy is a trademark of GE Healthcare
Pacific Blue and Pacific Orange are trademarks of Thermo Fisher Scientific