

LabNed PAP-Pen

Barrier Pen for IHC/IF

Ver. 06-2021

The LabNed PAP-Pen is a hydrophobic barrier pen that is detergent resistant and fixative resistant. It is widely used for IHC/IF applications for both FFPE and Frozen Tissue sections.

Advantages:

- No interference with fluorophores and antibodies
- Dries faster than similar products
- Can be washed off with Xylene
- Adheres to both Silane coated and plain slides
- Insensitive to most common detergents (Tween20, Triton X100 and X114)

Please check our website for ordering details.


* BarrierPaste only available on request



LabNed PAP-Pen Classic



Left to right: LabNed PAP-Pen Flat (7ml), Classic (7ml) and Classic (3ml)



Left to right: Line thickness of BarrierPaste*, PAP-Pen Flat and PAP-Pen Classic in all available colors

LabNed.com
exploring possibilities

