

# Puritan diagnostic swabs – made for the way you work.

**ImTec**  
DIAGNOSTICS NV

FOAM SWAB/  
TRANSPORT TUBE  
25-1506 1PF TT MC



1

Reclosable top  
for transport to  
the lab after  
sample collection

2

Molded base holds  
tip in center of tube  
to help avoid loss  
of sample material

3

Sterile foam swab  
and dry transport  
tube in one  
convenient unit

## WHY PURITAN DIAGNOSTIC SWABS?

We're known in the diagnostic field for **consistent quality – repeatable results – and dependable product availability**. Count on our swabs to deliver exactly what you need out in the field and in the lab.

As always, Puritan responds to market need with high-quality products. We are tuned into the diagnostic/healthcare world. We work hand-in-hand with our customers. By helping you solve challenges, we stay on the leading edge of our industry – creating the best single-use products to assist in reliable specimen collection and accurate diagnoses.


Our PurFlock swabs are a great example. The flock fibers are fully integrated with the handle by our proprietary manufacturing process so there are no heavy adhesives to interfere with bonding or with sample collection. The CDC has recommended sampling be done with synthetic tip materials and our flock swabs are ideal for those applications.

Take a look inside. We offer a large variety of swabs to meet just about any diagnostic need. And, please remember, custom orders are always welcome.



## HYDRAFLOCK® AND PURFLOCK® ULTRA FLOCKED SWABS

These patented specimen collection devices produce superior sample collection, rapid elution, and increased test sensitivity by combining a unique microstructure with innovative fibers proven to be more effective than older-style nylon products.



FEATURED SWABS: top to bottom # 3306-H / 3306-U, 3316-H / 3316-U, 3317-H / 3317-U, 3318-H / 3318-U

HydraFlock Item No.	PurFlock Ultra Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
3306-H	3306-U	Standard tip. Polystyrene handle with breakpoint at 80mm.	10/100, 3 lbs	152.654 x 2.54 mm tapered to 1.676 mm	17.526 x 5.08 mm 0.69 x 0.2"	153.416 mm 6.04"
25-3306-H	25-3306-U		10/50/1, 4 lbs	6.01 x 0.1" tapered to 0.066"		
3316-H	3316-U	Mini tip. Polystyrene handle available with breakpoint at 80mm or 100mm.	10/100, 3 lbs	152.4 x 2.489 mm tapered to 1.245 mm	18.415 x 3.429 mm 0.725 x 0.135"	153.416 mm 6.04"
25-3316-H	25-3316-U		10/50/1, 4 lbs	6 x 0.098" tapered to 0.049"		
3317-H	3317-U	Ultrafine tip. Polystyrene handle with breakpoint at 100mm.	10/100, 3 lbs	152.146 x 2.54 mm tapered to 0.914 mm	16.764 x 3.556 mm 0.66 x 0.14"	152.908 mm 6.02"
25-3317-H	25-3317-U		10/50/1, 4 lbs	5.99 x 0.1" tapered to 0.036"		
3318-H	3318-U	Micro ultrafine tip. Polystyrene handle with breakpoint at 100mm.	10/100, 3 lbs	152.146 x 2.54 mm tapered to 0.914 mm	7.925 x 3.556 mm 0.312 x 0.14"	153.162 mm 6.03"
25-3318-H	25-3318-U		10/50/1, 4 lbs	5.99 x 0.1" tapered to 0.036"		

## FOAM WOUND MEASUREMENT DEVICES

The soft foam-tipped applicator is flexible to conform to wound shape and features a metric rule for easy measuring. Perfect for forensics, emergency rooms, long-term care, and more. Choose sterile or non-sterile.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	1506-PFB DM	DM stick, polystyrene handle, foam tip.	20/50, 2 lbs	149.225 x 2.54 mm 5.875 x 0.1"	9.98 x 5.08 mm 0.393 x 0.2"	152.4 mm 6"
	25-1506 1PF DM		4/50/1, 2 lbs			

## FOAM-TIPPED APPLICATORS

These applicators are ideal for a variety of clinical and DNA applications. Tips are made of medical-grade foam. Choose from a range of shapes and handle materials.

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	non-sterile   sterile 25-1406 1PF 50	Foam tip and universal acetal handle. The small tip is ideal for patient comfort. Often used as a general swab for medical or veterinary diagnostics.	10/50/1, 4 lbs	152.298 x 2.489 mm tapered to 1.244 mm 5.996 x 0.098" tapered to 0.049"	13.894 x 3.581 mm 0.547 x 0.141"	152.4 mm 6"
	25-1506 1PF	Sterile foam swab for diagnostic or forensic sampling. Polystyrene handle.	10/100/1, 6 lbs	149.225 x 2.54 mm 5.875 x 0.1"	16.002 x 4.775 mm 0.63 x 0.188"	152.4 mm 6"
 product shown at 80%	25-1506 1PF TT	Swab 25-1506 1PF housed in a reclosable dry transport tube. Ideal for cell collection for rapid diagnostics and transport forensics.	24/30/1, 8 lbs	149.225 x 2.54 mm 5.875 x 0.1"	16.002 x 4.775 mm 0.63 x 0.188"	152.4 mm 6"
 product shown at 80%	25-1506 1PF TT MC	Swab 25-1506 1PF housed in a reclosable dry transport tube. Molded base holds the swab tip snugly in the tube center.	10/100/1, 10 lbs	149.225 x 2.54 mm 5.875 x 0.1"	16.002 x 4.775 mm 0.63 x 0.188"	152.4 mm 6"
	25-1605 1PSF RECT	Small rectangular foam tip has high collection capacity. For diagnostics, forensics, and DNA testing (including buccal swabs). Polypropylene handle.	10/50/1, 5 lbs	127 x 3.175 mm 5 x 0.125"	25.4 x 13.487 mm 1 x 0.531"	130.175 mm 5.125"
	25-1607 1PF SC	Foam tip, polypropylene handle. For detection of micro-organisms. Prominent thumb stop and large tip for maximum collection. Scored for easy detaching.	10/50/1, 5 lbs	127.635 x 2.997 mm 5.025 x 0.118"	25.4 x 11.125 mm 1 x 0.438"	130.81 mm 5.15"
	25-1805 1PF RND	Foam-tipped head offers large surface area for effective collection of small samples – ideal for luminol test of minute blood stains. Polypropylene handle.	10/50/1, 6 lbs	127 x 3.175 mm 5 x 0.125"	28.575 x 23.012 mm 1.125 x 0.906"	130.175 mm 5.125"
	2197	Cervisoft® has a non-abrasive foam tip for patient comfort. Ideal for sexual assault kits. Foam readily releases cells onto slides. Polystyrene handle.	10/100, 3 lbs	178.05 x 2.54 mm 7.01 x 0.1"	25.4 x 4.775 mm 1 x 0.188"	179.375 mm 7.062"
	25-2197		10/100/1, 7 lbs			

## RAYON-TIPPED APPLICATORS

Ideal for diagnostic specimen collection, wound treatment, and application of medicines. Tips are made from the finest USP grade of rayon fiber carded in our facility.

Choose an aluminum handle (shown here) or plastic, paper, or wood



FEATURED SWAB: # 25-801 R 50

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	800 R	Aluminum handle.	10/100, 1 lb	165.1 x 0.889 mm 6.5 x 0.035"	17.45 x 1.981 mm 0.687 x 0.078"	165.1 mm 6.547"
	25-800 R 50		10/50/1, 4 lbs	136.525 x 0.889 mm 5.375 x 0.035"		139.7 mm 5.5"
	25-801 R 50	Aluminum handle.	10/50/1, 4 lbs	136.525 x 0.635 mm 5.375 x 0.025"	17.45 x 1.981 mm 0.687 x 0.078"	139.7 mm 5.5"
	803-PR	Polystyrene handle.	10/10/100, 9 lbs	74.6 x 2.54 mm 2.937 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	76.2 mm 3"
	806-PR	Polystyrene handle.	10/10/100, 16 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	25-806 1PR		10/100/1, 6 lbs			
	25-806 1PAR	Fine paper handle.	10/100/1, 6 lbs	145.999 x 1.497 mm 5.748 x 0.059"	12.7 x 3.302 mm 0.5 x 0.13"	147.574 mm 5.81"
	806-WR	Wood handle.	10/10/100, 16 lbs	147.625 x 2.54 mm 5.812 x 0.1"	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"
	25-806 1WR		10/100/1, 6 lbs			
	25-806 2WR		10/100/2, 7 lbs			
 product shown at 60%	808	Paper handle; throat or OB/GYN swab.	10/50, 5 lbs	203.2 x 3.962 mm 8 x 0.156"	28.575 x 12.7 mm 1.125 x 0.5"	207.975 mm 8.188"
 product shown at 60%	25-808 1PR	Polystyrene handle; throat or OB/GYN swab.	10/50/1, 8 lbs	196.85 x 3.81 mm 7.75 x 0.15"	28.575 x 12.7 mm 1.125 x 0.5"	203.2 mm 8"
	25-808 2PR		10/50/2, 11 lbs			
 product shown at 60%	816	Polypropylene handle; proctoscopic applicator.	10/50, 8 lbs	403.225 x 3.81 mm 15.875 x 0.15"	28.575 x 12.7 mm 1.125 x 0.5"	406.4 mm 16"
 product shown at 60%	817	Polypropylene handle; proctoscopic applicator (extra binder for firmer tip).	10/50, 8 lbs	403.225 x 3.81 mm 15.875 x 0.15"	28.575 x 12.7 mm 1.125 x 0.5"	406.4 mm 16"
 product shown at 60%	818	Polystyrene handle, double-tipped, scored at four inches.	2,000 bulk, 11 lbs	196.85 x 3.81 mm 7.75 x 0.15"	28.575 x 11.099 mm 1.125 x 0.437"	203.2 mm 8"


## POLYESTER-TIPPED APPLICATORS

Synthetic applicator tips are made of medical-grade polyester fiber carded in our factory. They contain no fatty acids or latent sugars, making them ideal for specimen transport, lab work, and diagnostic kits.






	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	800 D	Aluminum handle.	10/100, 1 lb	165.1 x 0.889 mm 6.5 x 0.035"	17.45 x 1.981 mm 0.687 x 0.078"	166.294 mm 6.547"
	25-800 D 50		10/50/1, 4 lbs	136.525 x 0.889 mm 5.375 x 0.035"	17.45 x 1.981 mm 0.687 x 0.078"	139.7 mm 5.5"
	25-800 1PD 50	Ultrafine polystyrene handle.	10/50/1, 4 lbs	152.4 x 2.489 mm tapered to 1.32 mm 6 x 0.098" tapered to 0.052"	17.45 x 1.981 mm 0.687 x 0.078"	152.4 mm 6"
	800-PD ALUM	Aluminum/polystyrene handle.	10/100, 2 lbs	152.4 x 2.489 mm tapered to 0.889 mm 6 x 0.098" tapered to 0.035"	17.475 x 1.981 mm 0.688 x 0.078"	152.4 mm 6"
	25-800 1PD ALUM 50		10/50/1, 4 lbs			
	25-801 D 50	Aluminum handle.	10/50/1, 4 lbs	136.525 x 0.635 mm 5.375 x 0.025"	17.45 x 1.981 mm 0.687 x 0.078"	152.4 mm 6"
	803-PD	Polystyrene handle.	10/10/100, 10 lbs	74.599 x 2.54 mm 2.937 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	76.2 mm 3"
  	806-PD	Polystyrene handle.	10/10/100, 16 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	25-806 1PD		10/100/1, 6 lbs			
	25-806 2PD		10/100/2, 7 lbs			
	25-806 1PD TT	Polystyrene handle, reclosable dry transport tube.	24/30/1, 8 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	25-806 2PD TT		24/30/2, 9 lbs			
 product shown at 40%	25-806 1PP TT CAP	Polystyrene handle, reclosable dry transport tube; molded base holds swab tip snugly in center.	24/30/1, 8 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	806-WD	Wood handle.	10/10/100, 16 lbs	147.625 x 2.54 mm 5.812 x 0.1"	17.45 x 4.775 mm 0.687 x 0.188"	152.4 mm 6"
	25-806 1WD		10/100/1, 6 lbs			
	25-806 2WD		10/100/2, 7 lbs			
	25-806 1PD SOLID	Solid polystyrene handle.	10/100/1, 6 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	25-806 2PD SOLID		10/100/2, 8 lbs			
	25-826 2WD	Wood handle.	10/100/2, 6 lbs	147.625 x 2.108 mm 5.812 x 0.083"	12.7 x 3.048 mm 0.5 x 0.12"	152.4 mm 6"
	1000-PD	Polyphenylene Ether/polystyrene blend tapered handle, scored.	10/100, 2 lbs	150.114 x 2.159 mm tapered to 1.981 mm 5.91 x 0.085" tapered to 0.078"	17.45 x 1.981 mm 0.687 x 0.078"	152.4 mm 6"
	25-1000 1PD		10/100/1, 5 lbs			

## CALGISWAB® CALCIUM ALGINATE-TIPPED APPLICATORS

These specialty applicators are tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. Available with a variety of handle materials, including stainless steel upon request.



FEATURED SWAB: # 25-806 1PA



	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	non-sterile   sterile 25-800 A 50	Nasopharyngeal; aluminum handle.	10/50/1, 4 lbs	136.525 x 0.889 mm 5.375 x 0.035"	11.099 x 1.981 mm 0.437 x 0.078"	139.7 mm 5.5"
	25-801 A 50	Urethro-genital; aluminum handle.	10/50/1, 4 lbs	136.525 x 0.635 mm 5.375 x 0.025"	11.099 x 1.981 mm 0.437 x 0.078"	139.7 mm 5.5"
	806-PA 5M	Polystyrene handle.	5/10/100, 9 lbs	149.225 x 2.54 mm 5.875 x 0.1"	16.281 x 5.563 mm 0.641 x 0.219"	152.4 mm 6"
	25-806 1PA	Polystyrene handle.	10/100/1, 6 lbs	149.225 x 2.54 mm 5.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	152.4 mm 6"
	25-806 2PA		10/100/2, 7 lbs			
	806-WA 5M	Wood handle.	5/10/100, 9 lbs	147.625 x 2.54 mm 5.812 x 0.1"	16.281 x 5.563 mm 0.641 x 0.219"	152.4 mm 6"
	25-806 1WA		10/100/1, 6 lbs			
	25-806 2WA		10/100/2, 7 lbs			

## CAP-SHURE® COTTON-TIPPED APPLICATORS

Standard cotton-tipped applicator with aerated tip protector designed for convenient, reliable specimen, DNA or evidence collection. Comes as a unit with a reclosable, ventilated cap for air drying and transporting samples.



FEATURED SWAB: # 25-807 1PC EC

	Item No.	Item Description	Packaging pack, weight	Handle length x diameter	Tip length x diameter	Full Length
	non-sterile   sterile 25-806 1WC EC	Wood handle.	10/50/1, 7 lbs	147.625 x 2.54 mm 5.812 x 0.1"	15.875 x 4.775 mm 0.625 x 0.188"	152.4 mm 6"
	25-807 1PC EC	Polystyrene handle.	10/50/1, 9 lbs	174.625 x 2.54 mm 6.875 x 0.1"	17.45 x 5.156 mm 0.687 x 0.203"	177.8 mm 7"

## MAKE PURITAN SWABS YOUR OWN

Take advantage of our **custom options**, including handle color, printing, and variable break points.

### HANDLES ARE AVAILABLE IN YOUR CHOICE OF COLORS

The colors shown below are a few that customers have requested in the past. These are just examples – you can request any color you like. Colors make swab differentiation fast and easy – in the field, in the lab, and in diagnostic kits. You can also select a color to match your company or corporate palette.




### YOU CAN HAVE THE HANDLES IMPRINTED WITH YOUR LOGO OR TEXT

Customize any of our diagnostic swabs with a signature handle. Provide us with your logo, tag line, text, or other information and we'll imprint it right on the handle. For details about fonts and character counts, please get in touch with us.



Examples:



### ANOTHER OPTION FOR YOUR CONVENIENCE

Different applications require different handle lengths. Depending upon how you plan to use your Puritan diagnostic swab, you may have a specific break point in mind. Let us know and we'll build your break point into each molded handle.



### PACKAGING QUANTITY CODES

Our packaging codes work from the outside in. The first number represents the outermost package – this may be a case, a box, or other container.

Next is what you'll find when you open the outer packaging – bags, wrappers, tubes, etc.

Finally we get to the heart of the packaging – how many swabs are grouped together. Some examples follow.

10/50 = 10 bags of 50 swabs per bag, per case.

24/30/1 = 24 boxes of 30 swabs with each swab individually wrapped or in a tube.

10/100/1 = 10 boxes of 100 swabs with each swab individually wrapped or in a tube.

### HAVE QUESTIONS? CONTACT US.

We are always happy to talk about our products and help you choose those that are best for your applications.

Call us from the US and Canada at 800-321-2313, or from elsewhere at 207-876-3311.

If you like, send an email to [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com).





**Puritan**<sup>®</sup>  
DIAGNOSTICS



HydraFlock<sup>®</sup>



PurFlock<sup>®</sup>  
ULTRA



Puritan Medical  
A Green Company

Puritan is passionate about being green.

We think green and give energy efficiency and resource stewardship high priority. Below are just a few things that make us green. Visit our website for more.

- Puritan runs on **100% renewable energy** (wood process waste) for all manufacturing energy needs.
- We have an **ongoing energy upgrade initiative**.
- Puritan packaging materials are made up of at least **65% post-consumer recycled content**.

PURITAN MEDICAL PRODUCTS CO. LLC  
P.O. BOX 149, 31 SCHOOL STREET  
GUILFORD, MAINE 04443-0149

Tel 800-321-2313 (US and Canada)  
207-876-3311

Fax 800-323-4153 (US and Canada)  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

#### THE PURITAN ADVANTAGE™

Every one of our products in the Puritan and PurSwab lines comes with the Puritan Advantage™ built right in.

We are a US manufacturer, with R&D, fabrication, sterilization, packaging, shipping, QC and oversight, sales, and customer service under one roof. Our products are consistently reliable and readily available.

We are a global company. Our medical devices carry the CE mark.

ISO 9001:2008 ISO 13485:2003 CE

#### GET SOCIAL WITH PURITAN!



[facebook.com/PuritanUSA](https://facebook.com/PuritanUSA)



[twitter.com/PuritanUSA](https://twitter.com/PuritanUSA)



[linkedin.com/company/PuritanUSA](https://linkedin.com/company/PuritanUSA)



[youtube.com/PuritanUSA](https://youtube.com/PuritanUSA)